



Yeastern Biotech Co., Ltd

Product Use Limitation & Warranty

This product is intended to be used for life science research only. It has not been approved for drug or diagnostic purpose. YEASTERN's products should not be resold, modified for resale, or used to manufacture commercial products without written approval by YEASTERN. YEASTERN guarantees the performance of all products in the manner described in our protocol. The purchaser must determine the suitability of the product for its particular use. Should any product fail to perform satisfactorily due to any reason other than misuse, YEASTERN will replace it free of charge.

Ver. L0618

No part of these protocols may be reproduced in any form or by any mean, transmitted, or translated into a machine language without the permission of YEASTERN BIOTECH CO., LTD.

Address: 6F-3, 23 Lane 169, Kang Ning St.,
Shijr, Taipei, 22180 Taiwan.

Tel: +886-2-2695-3922 **Fax:** +886-2-2695-3979

Email: yeastern@yeastern.com.tw



YB Rapid Ligation Kit

Copyright© 2011 All rights reserved. Yeastern Biotech Co., Ltd.

Copyright© 2011 All rights reserved. Yeastern Biotech Co., Ltd.

Cat. No.
FYC003-100R

YB Rapid Ligation Kit

Concentration: 3 U/ μ l

Storage: -20 °C

Description

Yeastern's yT4 DNA ligase catalyzes the joining of two strands of DNA between the 5'-phosphate and the 3'-hydroxyl groups of adjacent nucleotides in either a cohesive-ended or blunt-ended termini. The enzyme also repairs single-strand nicks in duplex DNA, RNA or DNA/RNA hybrids. YB Rapid Ligation Kit is designed for efficient ligation of DNA inserts into vectors in just 5 minutes.

Component	Concentration	Volume
yT4 DNA ligase	3 U/ μ l	100 μ l
10X Ligation Buffer A		200 μ l
10X Ligation Buffer B		200 μ l

yT4 DNA ligase Storage Buffer:

20 mM Tris-HCl (pH 7.5), 50 mM KCl, 1 mM DTT, 0.1 mM EDTA and 50% glycerol.

10X Ligase Buffer A:

0.4 M Tris-HCl, 0.1 mM MgCl₂, 0.1 M DTT and 5 mM ATP (pH 5.0 at 25°C).

10X Ligase Buffer B:

A buffer contains an enhancer which dramatically increases ligation efficiency for blunt end DNA.

Unit Definition:

One unit of enzyme catalyzes the conversion of 1 nanomole of [³²P]Pi into Norit-adsorbable form in 20 min at 37°C (Weiss unit).

Standard Applications

- Joining double-stranded DNA with cohesive or blunt termini.
- Joining oligonucleotide linkers to blunt-ended DNA.
- Repairing nicks in duplex DNA, RNA or DNA-RNA hybrids.

Procedure

1. In a microcentrifuge tube prepare 5-10 μ l mix in ddH₂O or TE buffer if digested vector DNA and foreign DNA to be inserted.

2. Add the following components to the same tube:



	Cohesive Ends	Blunt Ends
vector: insert molar ratio	1 : 3	1 : 6
Vector fragments end conc.	3-30 fmol	15-60 fmol
Insert fragments end conc.	9-90 fmol	45-180 fmol
10X Ligation Buffer A	2 μ l	2 μ l
10X Ligation Buffer B	2 μ l	2 μ l
yT4 DNA ligase	1 μ l	1 μ l
ddH ₂ O to final volume of	20 μ l	20 μ l

3. Vortex the tube and spin down in microcentrifuge for 3-5 sec.
4. Incubate the mixture for 5-20 min. at 22°C.
5. If the insert fragment size > 3 kb, incubate the ligation mixture at 4°C overnight.
6. Use the mixture for transformation.

Note

1. yT4 DNA ligase is strongly inhibited by NaCl or KCl if the concentration exceeds 0.2 M.
2. 10X Ligation Buffer B greatly increase the rate of ligation of blunt-ended DNA.
3. Use equal or higher (up to 3-fold) molar concentration of insert DNA over vector DNA.
4. If the yield of ligation product is insufficient, prolong the reaction time (overnight).
5. Ligation reactions performed at lower temperatures require longer incubation time.
6. The performance of Ligation buffer A depends on the integrity of the ATP. Store the buffer in small aliquots at -20°C to minimize degradation of the ATP and DTT.
7. The DTT in the 10X Ligation Buffer A may precipitate upon freezing. If this occurs, vortex the buffer until the precipitate is completely dissolved in solution.
8. yT4 DNA Ligase is unstable on ice for a long period of time, please do not put it out of -20°C longer than 5-10 min.