

ECOS™ Sonic Competent Cells

I. Description

5/20 vials

Store at **-80°C**

ECOS™ competent cells is the fast transformation competent cells developed by Yeastern Biotech technology. Without additional gene modification, the process takes less than 6 minutes to yield the high efficiency without the need of SOC recovery (with ampicillin).

ECOS™ Sonic competent cells is the derivative of BL21(DE3) strain. It inherits the protein expression and fast-growing feature from its parent strain. The further deletion of *endA* and *recA* improve the quality of plasmid DNA and enable the plasmid extraction in 3-6 hours growing.


T&A Expression Kit (FYC201-10P) is recommended for ECOS™ Sonic; it combines the primary cloning of PCR product, and the subcloning for protein expression make the experiment easier and faster.

II. Genotype

F⁻ompT hsdS_B (r_B⁻m_B⁻) gal dcm (DE3)ΔendA ΔrecA

III. Efficiency

≥ 1*10⁸ cfu/μg



The concentration of plasmid purification will be largest after 3~6 hours incubation.



- PCR
- T&A Cloning
(Suggest: FYC201-10P)
or Other Cloning
- Transformation



- Plasmid
Identification



- Protein
Expression

IV. Transformation Procedure

We suggest you follow the 6-minute transformation protocol.

Important Note:

Ampicillin (Ap)	20 µg/ml
Kanamycin (Km)	25 µg/ml
Tetracycline (Tc)	25 µg/ml
Chloramphenicol (Cm)	20 µg/ml

1. The recommended concentration is based on freshly prepared antibiotics.
2. The use of expired antibiotics will result in insufficient medicinal properties and it is easy to produce pseudo-transformed colonies.
3. If an excessive concentration of antibiotics is used or a combination of multiple antibiotics is used for screening, the efficiency of the transformation will be significantly reduced.

